

Characteristic Groups in Rubber Extender and Processing Oils and Other Petroleum-Derived Oils by the Clay-Gel Absorption Chromatographic Method SAR-2007-01



Flexible & Solid Built

Method By Clay-Gel Absorption chromatographic

Test Method : ASTM D2007

Simplicity and safe design

SAR-2007-01 is built for with a strong base support to ensure the glassware are hold firmly due to the tall height. The clamps for holding the glassware can be configured to position at any preferred level and height by a convenient holder and screws. It features a complete set of glassware with clamps and heating mantle.

Technical Specifications

Unit Dimension : L1000 X B435 X H1500mm

Unit Weight 25 kg

Heating mantle power supply ; 230V50Hz

Set of complete glassware.

Optional recirculating chiller 230V 50Hz.

